

**INSTRUMEX B.O.D. INCUBATOR**

**Specifications:**

- Doubled walled inner chamber made of S.S and Outer chamber made of Mild Steel finished in powder coating epoxy paint.
- Full length inner acrylic door and outside metal door with magnetic gasket and lock
- Two removable S.S Tray at fix distance Unit consists of digital display temp. indicator cum controller with PT- 100 sensor
- Automatic chamber light and forced air circulation fan for uniform temp. Hermetically sealed compressor
- Castor wheel for easy move ability
- Temp. Range: 5 Deg. To 50 Deg. C accuracy +/-0.5 Deg. C.
- PID Controller Available On Request.

**Available in following sizes:**

**\*Also Special Sizes available on request.**

Capacity	Inner Dimension
3 Cu.ft. (90 Ltrs.)	45 * 45 * 50 Cms.
6 Cu.ft. (170 Ltrs.)	50 * 50 * 70 Cms.
10 Cu.ft. (280 Ltrs.)	55 * 55 * 90 Cms.
12 Cu.ft. (330 Ltrs.)	60 * 60 * 90 Cms.

